

# EOSINOPHIL DILUTING FLUID

**BIO LAB**  
**DIAGNOSTICS**  
ISO 9001:2015  
ISO 13485:2016  
CE

## PRINCIPLE

Eosinophils are selectively colored to pinkish orange while other white blood cells are uncolored. Diluted samples in eosinophil diluting fluid are charged to a Neubauer counting chamber and the pinkish colored cells are counted and calculated appropriately.

## REAGENTS COMPOSITION

Eosin Sodium	1 Gram/L
Acetic Acid	28.6mL/L
Distilled water	qs

## Working Reagent Preparation

Reagent is ready to use. May be filtered through watmann No 1 if precipitates or dust particles visible under microscope.

## STORAGE AND STABILITY.

Tight Capped Reagent is stable at room temperature until expiry date stated on the label.

## MATERIALS REQUIRED BUT NOT SUPPLIED.

1. Microscope with 45 x objectives
2. Neubauer counting chamber with cover glass.
3. 20  $\mu$ L pipette with tips.
4. 200  $\mu$ L pipette with tips.
5. Test tubes (75  $\times$  12 mm)
6. Petri dish with moist filter paper.


## SAMPLE COLLECTION AND HANDLING.

Whole blood collected in K<sub>3</sub> EDTA. Samples can be stored in 2-8°C for a week if the test need to perform at a later time. DO not Freez the sample below 2°C.

## PROCEDURE

1. Label test tubes with Sample ID
2. Pipette 200  $\mu$ L Eosinophil diluting fluid.
3. Pipette 20  $\mu$ L well mixed whole Blood sample.
4. Mix well and Charge in a Neubauer counting chamber.
5. Keep the charged chamber for 5 minutes under a petridish with moist filter paper.

## Neubauer counting chamber



6. Count only distinctly pinkish orange colored Leucocytes in all 4 corner WBC column squares under 45 x microscope.
7. Calculate absolute Eosinophil count as per result calculation column.

## RESULT CALCULATION

$$\text{Absolute Eosinophils / } \mu\text{L} = \frac{\text{TEC}}{4} \times \frac{220}{20} \times 10$$

$$\text{Absolute Eosinophils / } \mu\text{L} = \text{TEC} \times 27.5$$

Where TEC = Total Eosinophil Counted in 4 mm<sup>2</sup>  
 $\frac{220}{20}$  = Dilution of sample  
10 = Depth of counting chamber.

## EXPECTED VALUES

Absolute Eosinophils < 650 /  $\mu$ L  
Each laboratory should establish its own normal values.

## BIBLIOGRAPHY

1. Rider RF., Semin. Hematology 11.423 1974
2. John D Bauer., Numerical Evaluation of formed elements of blood., Gradwohl's clinical laboratory methods and diagnosis vol I page 797/9.

## BIOLAB DIAGNOSTICS (I) PVT. LTD.

J-245, MIDC, Tarapur, Boisar – 401 501, MS.  
E-mail : biolab@vsnl.com / www.biolabdiagnostics.com  
Customer Care : (+ 9122) 28088243