## CLINICAL SIGNIFICANCE

Sickle cells are RBC having abnormal hemoglobin called Hb-S which reduces the oxygen carrying capacity of Hemoglobin in blood circulation lead to various diseases called sickle cell anemia. The homozygous variants are with symptomatic where as heterozygous form of sickle cell diseases is asymptomatic called sickle cell trait. Patient with homozygous variants may have early mortality with vascular occlusion of multiple organs, severe hemolytic anemia and hypoxia. Patient with asymptomatic sickle cell disease may suffer sickle cell crisis such as hypoxia on anesthesia, during poorly pressurized air plain traveling, high altitude mountaineering or during pneumonia.

## PRINCIPLE

Hb-S precipitates in higher molar buffered medium in presence of Dithionate and detergent leaving non-sickling hemoglobin soluble which is free from turbidity. Whereas Dithionate reduced HbS get precipitated and form translucent in the test system.

## Kit Content

1. Reagent R1	100 mL
2. Reagent R2	1 bottle
3. Test Tubes and caps	50 nos
4. Reagent transfer pipette	1 No
5. Rapid Sickle Cell Reader rack	1 No
6. Pack insert	1 No

## REAGENTS COMPOSITION

### Reagent 1 : R1

di-Potassium phosphate	2.3 mol/L
Detergent	0.4 g/L
pH7.1(±0.1) Contain Stabilizers and preservatives.	

### Reagent 2 : R2

Dithionate	57 mmol/L
White Saponins	1 %

## Working Reagent Preparation

Transfer content of one Vial of R2 to a bottle of R1 and mix gently for 15 mins. Make sure that entire content of R2 is transferred to R1. Label date of reconstitution.

## STORAGE AND STABILITY

When stored at room temperature and protected from light, the reagent is stable 12 months from the manufacturing date stated on the label. Working reagent is stable for 15 days at ambient room temperature, 6 months at 2-8°C. **Do not freeze or expose to light.** Discard deteriorated working reagent that turns yellowish when mixed with blood.

## SAMPLES

No Prior patient preparation is needed. Use fresh blood collected from clean finger prick or in EDTA. More accuracy may be achieved if the test is performed with Packed cells by leaving the sample for 30 minutes undisturbed; RBC sediment can be collected from the bottom of sample container. Sample can be stored for a week at 4 to 8 °C if the test needs to perform at a later time. (All samples should be handled as potential infective agents as no laboratory methods make conclusive finding for its safety. Therefore, adequate protective laboratory measures should be taken while handling such materials).


**Note.** Strong lipemic or hemolysed blood sample should not be used.

## PROCEDURE

1. Label each test tube with patient ID and pour 2 mL Working reagent in all the tubes using reagent transfer pipette.
2. Collect one drop (20 µL) of patient blood using a pipette tip connected to an appropriate adaptor or pipette. Mix well and wait for 10 minutes at room temperature ( stable for an hour), mix and interpret the result as given under.

**\* Sample having very low (< 7g%) hemoglobin perform test with 2 drops (40 µL) in place of one drop (20 µL) of blood.**

## INTERPRETATION OF RESULT

Mix the tubes and observe a Negative test as transparent to see the background in Rapid Sickle Cell Reader rack. Where as the background is not visible in positive sample (HbAS or SS) given as above.

## RESULTS VISUAL

Observe the result after 10 minutes of sample mixing. A positive sickling test indicates relatively less transparent than negative control while viewing in Rapid Sickle Cell Reader rack (refer interpretation of result column). Run a known negative & positive blood sample of similar hemoglobin concentration as negative and positive control for comparison. A positive or doubtful result must be reconfirmed with Hb Electrophoresis.

## INTERPRETIVE VALUES

### Strong Turbid

Homozygous for Sickling Hb - > 75 % of Total Hb %

### Weak Turbid

Heterozygous (Trait) for Sickling Hb - > 40 to < 70 % Total Hb %

### Faint Turbid

Insignificant - < 25 % Total Hb %  
(Re-assay with 2 drops blood)

**Note:** It is recommended for each laboratory to establish and maintain its own reference values. The above values are only an indication and to be confirmed through Electrophoresis.

## LIMITATION

Solubility test is only a screening method. This method is not a substitute to Hb Electrophoresis. A variant of Hb C – Harlem interfere this test. Recent blood transfusion may alter the results. Higher Hb F lead to false positive test result and very anemic or sample with less RBC may lead to false negative observation. Borderline Hb AS or Hb SS need to be interpreted with professional experience. Proper visibility needs appropriate illumination. Interpretation of result at inadequate light may lead to false results. Sample having very low (< 7g%) hemoglobin may lead to false negative. Repeat test using 2 drops (40 µL) in place of one drop (20 µL) of blood.

## WARNING

This reagent system is for in vitro use Only. This reagent system contains preservatives and components that have not established for safety if contacted on broken skin or eye or taken orally. In case of such incidents wash off with plenty of water and consult a physician.

## QUALITY CONTROL

To ensure adequate quality control, each analysis should be tested against known positive and negative blood samples.

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