

# BIOLYSE-3™

## Laboratory Cleaning Solution



ISO 9001:2015  
ISO 13485:2016  
CE

**BIOLYSE-3™** is a liquid disinfectant, Deodorant, Detergent with broad spectrum bactericidal, Fungicidal and viricidal actions effective against bacteria, fungi and viruses. **BIOLYSE-3™** is an ideal disinfectant protecting yourself and simultaneously your patients.

### APPLICATION

As **BioLyse-3™** is containing surfactant and Meta Cresol. It could be used in all places where cleaning is required such as floors, benches, sanitary items where infectious contaminants are susceptible.

### SPECIFICATION

pH - 7 (Neutral)  
Concentration 10x  
Soluble in water.

### MODE OF USE

Make 10% solution by diluting 100 mL **BIOLYSE-3™** into a liter of water. Use working solution for all cleaning purpose.

### BIOLYSE-3™ ENSURES

Aseptic Environment.  
Protection From Microbes.  
Cleaning and Deodorizing.

### BIOCLEAN ENSURES

Best cleaning on any surface. It is a stable, efficient surface active agent with a narrow range of pH variation. A 5% solution has a pH 6.5 to 7.5 that is near neutral pH.

Non-Corrosive, Non-toxic

### PRECAUTIONS

**BIOCLEAN** is for laboratory use only.  
Not for food, drug or household use.  
Avoid contact with skin, eyes.  
Keep away from children.

### PACKING

1000 mL Plastic bottles and  
5000 mL Plastic Jerry cans.

### BIO LAB DIAGNOSTICS (I) PVT. LTD.

J-245, MIDC, Tarapur, Boisar – 401 501, MS.

E-mail : [biolab@vsnl.com](mailto:biolab@vsnl.com)

[www.biolabdiagnostics.com](http://www.biolabdiagnostics.com)

Customer Care : (+ 9122) 28088243